

A decorative graphic in the top-left corner consisting of several orange circles of varying sizes, some solid and some with a white center, arranged in a cluster.

## Tobii Pro Eye Tracker Manager Release Notes

### Version 1.6

**Release date: 11 April 2018**

- Support for choosing "human" or "macaque" mode with eye tracker Tobii Pro Spectrum
- Bug fixes

### Version 1.5

**Release date: 20 February 2018**

- Support for Windows 7
- Support for Windows 32 bits
- Added release notes for Eye Tracker Manager application updates
- Added notification for Calibration and Display setup changes
- UX and UI improvements
- Possibility to include eye images in diagnostic file on specific eye trackers
- Bug fix for off-centered icons in Mac

### Version 1.4.1

**Release date: 11 December 2017**

- Fixed a bug related to Display setup mapping calculations on Tobii Pro TX300 and Tobii Pro Spectrum

### Version 1.4

**Release date: 22 November 2017**

- Improved user calibration
- Possibility to download and install a new available firmware of Tobii Pro Spectrum by giving notification
- Possibility to download and install a new version of the Eye Tracker Manager Application by giving notification
- Support for Mobile Device Stand display setup
- Chinese and Japanese language options
- Support for detecting Tobii Pro VR with information about serial number and firmware

## Version 1.3

**Release date: 11 October 2017**

- Improved user interface
- Eye Images visualization for Tobii Pro Spectrum
- Support for storing Tobii Pro Upgrade Key on Tobii Eye Tracker 4C
- Support for reusing some of the features of Tobii Pro Eye Tracker Manager in your own applications for SDK users (More information: <http://developer.tobii.com/eyetrackermanager/etm-sdk-integration.html>)

## Version 1.2

**Release date: 5 July 2017**

- Eye Tracker Manager installer is available for Linux
- Improved user calibration
- Support for saving Eye Tracker's diagnostic files
- Support for standalone display setup of Tobii Pro TX300, Tobii Pro X2 and Tobii Pro X3

## Version 1.1

**Release date: 24 May 2017**

- Improved user calibration
- Support for display setup of Tobii Pro X2 and Tobii Pro X3 and Tobii Eye Tracker 4C
- Track status visualization of the user's eyes
- Fixed a bug where the application Window in some corners got stuck in full screen

## Version 1.0

**Release date: 12 May 2017**

- Support for updating Tobii Pro spectrum firmware
- Support for changing eye tracking frequency
- Support for display setup of Tobii Pro Spectrum
- Support for detecting Tobii Pro Eye Trackers with information about serial number and firmware

